





LightWorks LW2

Troubleshooting





LightWorks LW2

Troubleshooting

Cannot switch the operational side of the Hand Paddle.

Ensure both ends of the LightWorks Cord are securely plugged into the unit and the Hand Paddle.

Yes

Are the Hand Paddle lights on?

No

Yes

Can the side in use be changed by pushing the Hand Paddle button?

No

The unit appears to be working properly.

There may be a fault with the unit. Contact Customer Service at 1.800.224.0242 or repairs@sota.com

Cannot switch modes.

There may be a fault with the unit. Contact Customer Service at 1.800.224.0242 or repairs@sota.com

When turned on, the Hand Paddle lights immediately flash and the LightWorks unit beeps and turns off.

Indicates the Wall Adaptor in use does not meet the minimum requirements to power the LightWorks.

The Wall Adaptor should be 12 Volts DC with a Center Positive Pin that supplies at least 1.5 Amps.